

# CAMK2N1 Antibody, Biotin conjugated



PACO63630

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Potent and specific inhibitor of CaM-kinase II (CAMK2).

**Gene ID:**

CAMK2N1

**Uniprot**

Q7Z7J9

**Synonyms:**

Calcium/calmodulin-dependent protein kinase II inhibitor 1 (CaMKII inhibitory protein alpha) (CaMKIIN-alpha), CAMK2N1

**Immunogen:**

Peptide sequence from Human Calcium/calmodulin-dependent protein kinase II inhibitor 1 protein (9-27AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A