

CAMK2N2 Antibody, Biotin conjugated



PACO63634

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Potent and specific cellular inhibitor of CaM-kinase II (CAMK2). Traps Ca(2+)/calmodulin on CAMK2. May play an important role in the regulation of cell growth when overexpressed in colon adenocarcinoma LoVo cells. Traps Ca(2+)/calmodulin on CAMK2.

Gene ID:

CAMK2N2

Uniprot

Q96S95

Synonyms:

Calcium/calmodulin-dependent protein kinase II inhibitor 2 (CaM-KII inhibitory protein) (CaM-KIIN), CAMK2N2

Immunogen:

Peptide sequence from Human Calcium/calmodulin-dependent protein kinase II inhibitor 2 protein (10-28AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A