

# CLCN5 Antibody, FITC conjugated



PACO63641

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Proton-coupled chloride transporter. Functions as antiport system and exchanges chloride ions against protons. Important for normal acid, fication of the endosome lumen. May play an important role in renal tubular function.

**Gene ID:**

CLCN5

**Uniprot**

P51795

**Synonyms:**

H(+)/Cl(-) exchange transporter 5 (Chloride channel protein 5) (ClC-5) (Chloride transporter ClC-5), CLCN5, CLCK2

**Immunogen:**

Recombinant Human H(+)/Cl(-) exchange transporter 5 protein (553-746AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A