

CLCN5 Antibody, Biotin conjugated



PACO63642

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Proton-coupled chloride transporter. Functions as antiport system and exchanges chloride ions against protons. Important for normal acid, fication of the endosome lumen. May play an important role in renal tubular function.

Gene ID:

CLCN5

Uniprot

P51795

Synonyms:

H(+)/Cl(-) exchange transporter 5 (Chloride channel protein 5) (ClC-5) (Chloride transporter ClC-5), CLCN5, CLCK2

Immunogen:

Recombinant Human H(+)/Cl(-) exchange transporter 5 protein (553-746AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A