

# COL6A1 Antibody, FITC conjugated



PACO63645

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Collagen VI acts as a cell-binding protein.

**Gene ID:**

COL6A1

**Uniprot**

P12109

**Synonyms:**

Collagen alpha-1(VI) chain, COL6A1

**Immunogen:**

Peptide sequence from Human Collagen alpha-1(VI) chain protein (64-80AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A