

# GNAT3 Antibody, HRP conjugated



PACO63656

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

GNAT3

**Uniprot**

A8MTJ3

**Synonyms:**

Guanine nucleotide-binding protein G(t) subunit alpha-3 (Gustducin alpha-3 chain), GNAT3

**Immunogen:**

Recombinant Human Guanine nucleotide-binding protein G(t) subunit alpha-3 protein (2-354AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A