

GNAT3 Antibody, FITC conjugated



PACO63657

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

GNAT3

Uniprot

A8MTJ3

Synonyms:

Guanine nucleotide-binding protein G(t) subunit alpha-3 (Gustducin alpha-3 chain), GNAT3

Immunogen:

Recombinant Human Guanine nucleotide-binding protein G(t) subunit alpha-3 protein (2-354AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A