

PITPNA Antibody, Biotin conjugated



PACO63694

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the transfer of PtdIns and phosphatidylcholine between membranes.

Gene ID:

PITPNA

Uniprot

Q00169

Synonyms:

Phosphatidylinositol transfer protein alpha isoform (PI-TP-alpha) (PtdIns transfer protein alpha) (PtdInsTP alpha), PITPNA, PITPN

Immunogen:

Recombinant Human Phosphatidylinositol transfer protein α isoform protein (1-270AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A