

S100A7 Antibody, HRP conjugated



PACO63708

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

azuophil granule lumen, cytoplasm, cytosol, endoplasmic reticulum, extracellular exosome, extracellular region, extracellular space, focal adhesion, nucleus, calcium ion binding.

Gene ID:

S100A7

Uniprot

P31151

Synonyms:

Protein S100-A7 (Psoriasin) (S100 calcium-binding protein A7), S100A7, PSOR1 S100A7C

Immunogen:

Recombinant Human Protein S100-A7 protein (2-101AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A