

# S100A7 Antibody, FITC conjugated



PACO63709

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

azurophil granule lumen, cytoplasm, cytosol, endoplasmic reticulum, extracellular exosome, extracellular region, extracellular space, focal adhesion, nucleus, calcium ion binding.

**Gene ID:**

S100A7

**Uniprot**

P31151

**Synonyms:**

Protein S100-A7 (Psoriasin) (S100 calcium-binding protein A7), S100A7, PSOR1 S100A7C

**Immunogen:**

Recombinant Human Protein S100-A7 protein (2-101AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A