

**Uncharacterized protein Antibody, HRP  
conjugated  
PACO63740**



---

**Product Information**

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

Uncharacterized protein

**Uniprot**

E7EVH7

**Synonyms:**

Uncharacterized protein

**Immunogen:**

Recombinant Human Uncharacterized protein (502-717AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A