

Uncharacterized protein Antibody, Biotin conjugated
PACO63742



Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uncharacterized protein

Uniprot

E7EVH7

Synonyms:

Uncharacterized protein

Immunogen:

Recombinant Human Uncharacterized protein (502-717AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A