

# TRAPPC8 Antibody, HRP conjugated



PACO63744

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May be involved in endoplasmic reticulum to Golgi apparatus trafficking at a very early stage.

**Gene ID:**

TRAPPC8

**Uniprot**

Q9Y2L5

**Synonyms:**

Trafficking protein particle complex subunit 8 (Protein TRS85 homolog), TRAPPC8, KIAA1012

**Immunogen:**

Recombinant Human Trafficking protein particle complex subunit 8 protein (1205-1435AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A