

TRAPPC8 Antibody, Biotin conjugated



PACO63746

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in endoplasmic reticulum to Golgi apparatus trafficking at a very early stage.

Gene ID:

TRAPPC8

Uniprot

Q9Y2L5

Synonyms:

Trafficking protein particle complex subunit 8 (Protein TRS85 homolog), TRAPPC8, KIAA1012

Immunogen:

Recombinant Human Trafficking protein particle complex subunit 8 protein (1205-1435AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A