

CDK17 Antibody, Biotin conjugated



PACO63750

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in terminally differentiated neurons. Has a Ser/Thr-phosphorylating activity for histone H1.

Gene ID:

CDK17

Uniprot

Q00537

Synonyms:

Cyclin-dependent kinase 17 (EC 2.7.11.22) (Cell division protein kinase 17) (PCTAIRE-motif protein kinase 2) (Serine/threonine-protein kinase PCTAIRE-2), CDK17, PCTAIRE2, PCK2

Immunogen:

Recombinant Human Cyclin-dependent kinase 17 protein (337-523AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A