

# N Antibody, HRP conjugated



PACO63756

---

## Product Information

**Size:**

50ul

**Reactivity:**

Bovine respiratory syncytial virus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Encapsidates the genome, protecting it from nucleases. The nucleocapsid (NC) has a helical structure. The encapsidated genomic RNA is termed the NC and serves as template for transcription and replication. During replication, encapsidation by N is coupled to RNA synthesis and all replicative products are resistant to nucleases.

**Gene ID:**

N

**Uniprot**

P35943

**Synonyms:**

Nucleoprotein (Protein N) (Nucleocapsid protein), N

**Immunogen:**

Recombinant Bovine respiratory syncytial virus Nucleoprotein protein (1-391AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A