

PCOLCE Antibody, Biotin conjugated



PACO63794

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds to the C-terminal propeptide of type I procollagen and enhances procollagen C-proteinase activity.

Gene ID:

PCOLCE

Uniprot

Q15113

Synonyms:

Procollagen C-endopeptidase enhancer 1 (Procollagen COOH-terminal proteinase enhancer 1) (PCPE-1) (Procollagen C-proteinase enhancer 1) (Type 1 procollagen C-proteinase enhancer protein) (Type I procollagen COOH-terminal proteinase enhancer), PCOLCE, PCPE1

Immunogen:

Recombinant Human Procollagen C-endopeptidase enhancer 1 protein (169-449AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A