

GRM3 Antibody, Biotin conjugated



PACO63798

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

G-protein coupled receptor for glutamate. Ligand binding causes a conformation change that triggers signaling via guanine nucleotide-binding proteins (G proteins) and modulates the activity of down-stream effectors. Signaling inhibits adenylate cyclase activity.

Gene ID:

GRM3

Uniprot

Q14832

Synonyms:

Metabotropic glutamate receptor 3 (mGluR3), GRM3, GPRC1C MGLUR3

Immunogen:

Recombinant Human Metabotropic glutamate receptor 3 protein (46-354AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A