

# omp Antibody



PACO63803

---

## Product Information

**Size:**

50ul

**Reactivity:**

Rickettsia japonica

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:1000-1:5000

**Protein Background:****Gene ID:**

omp

**Uniprot**

Q52764

**Synonyms:**

17 kDa surface antigen, omp

**Immunogen:**


Recombinant Rickettsia japonica 17 kDa surface antigen protein (20-159AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Western Blot. Positive WB detected in Recombinant protein. All lanes: omp antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 32 kDa. Observed band size: 32 kDa.