

# DNER Antibody, Biotin conjugated



PACO63822

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Activator of the NOTCH1 pathway. May mediate neuron-glia interaction during astrocytogenesis.

**Gene ID:**

DNER

**Uniprot**

Q8NFT8

**Synonyms:**

Delta and Notch-like epidermal growth factor-related receptor, DNER, BET

**Immunogen:**

Peptide sequence from Human Delta and Notch-like epidermal growth factor-related receptor protein (150-168AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A