

RORB Antibody, Biotin conjugated



PACO63830

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

RORB

Uniprot

Q92753

Synonyms:

Nuclear receptor ROR-beta (Nuclear receptor RZR-beta) (Nuclear receptor subfamily 1 group F member 2) (Retinoid-related orphan receptor-beta), RORB, NR1F2 RZRB

Immunogen:

Recombinant Human Nuclear receptor ROR-β protein (1-269AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A