

# MTMR9 Antibody, HRP conjugated



PACO63848

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Probable pseudophosphatase. Contains a Gly residue instead of a conserved Cys residue in the dsPTPase catalytic loop which renders it catalytically inactive as a phosphatase (Potential).

**Gene ID:**

MTMR9

**Uniprot**

Q96QG7

**Synonyms:**

Myotubularin-related protein 9, MTMR9, C8orf9 MTMR8

**Immunogen:**

Recombinant Human Myotubularin-related protein 9 protein (322-549AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A