

STYX Antibody, Biotin conjugated



PACO63858

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Catalytically inactive phosphatase. Acts as a nuclear anchor for MAPK1/MAPK3 (ERK1/ERK2). Modulates cell-fate decisions and cell migration by spatiotemporal regulation of MAPK1/MAPK3 (ERK1/ERK2). Seems to play a role in spermiogenesis.

Gene ID:

STYX

Uniprot

Q8WUJ0

Synonyms:

Serine/threonine/tyrosine-interacting protein (Protein tyrosine phosphatase-like protein), STYX

Immunogen:

Recombinant Human Serine/threonine/tyrosine-interacting protein (1-223AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A