

CAPN10 Antibody, HRP conjugated



PACO63880

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Calcium-regulated non-lysosomal thiol-protease which catalyze limited proteolysis of substrates involved in cytoskeletal remodeling and signal transduction. May play a role in insulin-stimulated glucose uptake.

Gene ID:

CAPN10

Uniprot

Q9HC96

Synonyms:

Calpain-10 (EC 3.4.22) (Calcium-activated neutral proteinase 10) (CANP 10), CAPN10, KIAA1845

Immunogen:

Recombinant Human Calpain-10 protein (1-287AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A