

CLDND1 Antibody, Biotin conjugated



PACO63898

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cell surface.

Gene ID:

CLDND1

Uniprot

Q9NY35

Synonyms:

Claudin domain-containing protein 1 (Membrane protein GENX-3745), CLDND1, C3orf4

Immunogen:

Recombinant Human Claudin domain-containing protein 1 protein (26-140AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A