

CSNK1G3 Antibody, HRP conjugated



PACO63920

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Serine/threonine-protein kinase. Casein kinases are operationally defined by their preferential utilization of acid, c proteins such as caseins as substrates. It can phosphorylate a large number of proteins. Participates in Wnt signaling. Regulates fast synaptic transmission mediated by glutamate.

Gene ID:

CSNK1G3

Uniprot

Q9Y6M4

Synonyms:

Casein kinase I isoform gamma-3 (CKI-gamma 3) (EC 2.7.11.1), CSNK1G3

Immunogen:

Recombinant Human Casein kinase I isoform gamma-3 protein (1-204AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A