

CALB1 Antibody, Biotin conjugated



PACO63950

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Buffers cytosolic calcium. May stimulate a membrane Ca(2+)-ATPase and a 3',5'-cyclic nucleotide phosphodiesterase.

Gene ID:

CALB1

Uniprot

P05937

Synonyms:

Calbindin (Calbindin D28) (D-28K) (Vitamin D-dependent calcium-binding protein, avian-type), CALB1, CAB27

Immunogen:

Recombinant Human Calbindin protein (2-261AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A