

# OPCML Antibody, FITC conjugated



PACO63969

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Binds opioids in the presence of acid, c lipids; probably involved in cell contact.

**Gene ID:**

OPCML

**Uniprot**

Q14982

**Synonyms:**

Opioid-binding protein/cell adhesion molecule (OBCAM) (OPCML) (Opioid-binding cell adhesion molecule) (IgLON family member 1), OPCML, IGLON1 OBCAM

**Immunogen:**

Recombinant Human Opioid-binding protein/cell adhesion molecule protein (28-322AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A