

DNASE2B Antibody, Biotin conjugated



PACO63978

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Hydrolyzes DNA under acid, c conditions. Does not require divalent cations for activity. Participates in the degradation of nuclear DNA during lens cell differentiation.

Gene ID:

DNASE2B

Uniprot

Q8WZ79

Synonyms:

Deoxyribonuclease-2-beta (EC 3.1.22.1) (DNase II-like acid, DNase) (DNase2-like acid, DNase) (Deoxyribonuclease II beta) (DNase II beta) (Endonuclease DLAD), DNASE2B, DLAD

Immunogen:

Recombinant Human Deoxyribonuclease-2-β protein (68-281AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A