

DPEP1 Antibody, Biotin conjugated



PACO63982

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Hydrolyzes a wide range of dipeptides. Implicated in the renal metabolism of glutathione and its conjugates. Converts leukotriene D4 to leukotriene E4; it may play an important role in the regulation of leukotriene activity.

Gene ID:

DPEP1

Uniprot

P16444

Synonyms:

Dipeptidase 1 (EC 3.4.13.19) (Dehydropeptidase-I) (Microsomal dipeptidase) (Renal dipeptidase) (hRDP), DPEP1, MDP RDP

Immunogen:

Recombinant Human Dipeptidase 1 protein (23-252AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A