

EXOC3 Antibody, FITC conjugated



PACO63989

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Component of the exocyst complex involved in the docking of exocytic vesicles with fusion sites on the plasma membrane.

Gene ID:

EXOC3

Uniprot

O60645

Synonyms:

Exocyst complex component 3 (Exocyst complex component Sec6), EXOC3, SEC6, SEC6L1

Immunogen:

Recombinant Human Exocyst complex component 3 protein (410-756AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A