

IKBIP Antibody, Biotin conjugated



PACO64026

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Target of p53/TP53 with pro-apoptotic function.

Gene ID:

IKBIP

Uniprot

Q70UQ0

Synonyms:

Inhibitor of nuclear factor kappa-B kinase-interacting protein (I kappa-B kinase-interacting protein) (IKBKB-interacting protein) (IKK-interacting protein), IKBIP, IKIP

Immunogen:

Recombinant Human Plexin-B2 protein (121-350AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A