

# TFF3 Antibody, HRP conjugated



PACO64032

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Involved in the maintenance and repair of the intestinal mucosa. Promotes the mobility of epithelial cells in healing processes (motogen).

**Gene ID:**

TFF3

**Uniprot**

Q07654

**Synonyms:**

Trefoil factor 3 (Intestinal trefoil factor) (hITF) (Polypeptide P1. B) (hP1. B), TFF3, ITF TFI

**Immunogen:**

Recombinant Human Zinc finger protein 346 protein (22-80AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A