

TFF3 Antibody, Biotin conjugated



PACO64034

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in the maintenance and repair of the intestinal mucosa. Promotes the mobility of epithelial cells in healing processes (motogen).

Gene ID:

TFF3

Uniprot

Q07654

Synonyms:

Trefoil factor 3 (Intestinal trefoil factor) (hITF) (Polypeptide P1. B) (hP1. B), TFF3, ITF TFI

Immunogen:

Recombinant Human Zinc finger protein 346 protein (22-80AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A