

# UBB Antibody, Biotin conjugated



PACO64066

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

UBB

**Uniprot**

P0CG47

**Synonyms:**

Polyubiquitin-B [Cleaved into: Ubiquitin], UBB

**Immunogen:**

Recombinant Human Polyubiquitin-B protein (103-228AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A