

# CXCL12 Antibody, HRP conjugated



PACO64068

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CXCL12

**Uniprot**

P48061

**Synonyms:**

Stromal cell-derived factor 1 (SDF-1) (hSDF-1) (C-X-C motif chemokine 12) (Intercrine reduced in hepatomas) (IRH) (hIRH) (Pre-B cell growth-stimulating factor) (PBSF) [Cleaved into: SDF-1-beta(3-72), SDF-1-alpha(3-67)], CXCL12, SDF1 SDF1A SDF1B

**Immunogen:**

Recombinant Human Stromal cell-derived factor 1 protein (22-89AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A