

KCNK15 Antibody, HRP conjugated



PACO64120

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probable potassium channel subunit. No channel activity observed in heterologous systems. May need to associate with another protein to form a functional channel.

Gene ID:

KCNK15

Uniprot

Q9H427

Synonyms:

Potassium channel subfamily K member 15 (acid, sensitive potassium channel protein TASK-5) (TWIK-related acid, sensitive K(+) channel 5) (Two pore potassium channel KT3.3) (Two pore K(+) channel KT3.3), KCNK15, TASK5

Immunogen:

Recombinant Human Potassium channel subfamily K member 15 protein (244-330AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A