

APOM Antibody, HRP conjugated



PACO64144

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probably involved in lipid transport. Can bind sphingosine-1-phosphate, myristic acid, palmitic acid, and stearic acid, retinol, all-trans-retinoic acid, and 9-cis-retinoic acid, .

Gene ID:

APOM

Uniprot

O95445

Synonyms:

Apolipoprotein M (Apo-M) (ApoM) (Protein G3a), APOM, G3A NG20

Immunogen:

Recombinant Human Apolipoprotein M protein (1-188AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A