

# APOM Antibody, Biotin conjugated



PACO64146

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Probably involved in lipid transport. Can bind sphingosine-1-phosphate, myristic acid, palmitic acid, and stearic acid, retinol, all-trans-retinoic acid, and 9-cis-retinoic acid, .

**Gene ID:**

APOM

**Uniprot**

O95445

**Synonyms:**

Apolipoprotein M (Apo-M) (ApoM) (Protein G3a), APOM, G3A NG20

**Immunogen:**

Recombinant Human Apolipoprotein M protein (1-188AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A