

# CCDC74B Antibody



PACO64187

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:1000-1:5000

**Protein Background:**

**Gene ID:**

CCDC74B

**Uniprot**

Q96LY2

**Synonyms:**

Coiled-coil domain-containing protein 74B, CCDC74B

**Immunogen:**


Recombinant Human Coiled-coil domain-containing protein 74B protein (1-257AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Western Blot. Positive WB detected in: PC-3 whole cell lysate, HeLa whole cell lysate, MCF-7 whole cell lysate. All lanes: CCDC74B antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 42, 36 kDa. Observed band size: 42 kDa.