

CHCHD5 Antibody, FITC conjugated



PACO64197

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

mitochondrial intermembrane space, aerobic respiration.

Gene ID:

CHCHD5

Uniprot

Q9BSY4

Synonyms:

Coiled-coil-helix-coiled-coil-helix domain-containing protein 5, CHCHD5, C2orf9

Immunogen:

Recombinant Human Coiled-coil-helix-coiled-coil-helix domain-containing protein 5 protein (1-110AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A