

CLCN6 Antibody, HRP conjugated



PACO64200

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Chloride transport protein, initially identified as voltage-gated chloride channel. The presence of the conserved gating glutamate residues suggests that it functions as antiporter.

Gene ID:

CLCN6

Uniprot

P51797

Synonyms:

Chloride transport protein 6 (Chloride channel protein 6) (CIC-6), CLCN6, KIAA0046

Immunogen:

Recombinant Human Chloride transport protein 6 protein (572-869AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A