

# CLDN11 Antibody, FITC conjugated



PACO64205

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Plays a major role in tight junction-specific obliteration of the intercellular space, through calcium-independent cell-adhesion activity.

**Gene ID:**

CLDN11

**Uniprot**

O75508

**Synonyms:**

Claudin-11 (Oligodendrocyte-specific protein), CLDN11, OSP OTM

**Immunogen:**

Recombinant Human Claudin-11 protein (23-82AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A