

FAM24B Antibody, Biotin conjugated



PACO64230

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

FAM24B

Uniprot

Q8N5W8

Synonyms:

Protein FAM24B, FAM24B

Immunogen:

Recombinant Human Protein FAM24B protein (22-94AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A