

FBXO36 Antibody, Biotin conjugated



PACO64238

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex.

Gene ID:

FBXO36

Uniprot

Q8NEA4

Synonyms:

F-box only protein 36, FBXO36, FBX36

Immunogen:

Recombinant Human F-box only protein 36 protein (1-188AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A