

PACO64239

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:1000-1:5000

**Protein Background:**

Binds directly to activin and functions as an activin antagonist. Specific inhibitor of the biosynthesis and secretion of pituitary follicle stimulating hormone (FSH).

**Gene ID:**

FST

**Uniprot**

P19883

**Synonyms:**

Follistatin (FS) (Activin-binding protein), FST

**Immunogen:**


Recombinant Human Follistatin protein (30-344AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Western Blot. Positive WB detected in: U251 whole cell lysate. All lanes: FST antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 39, 35 kDa. Observed band size: 39 kDa.