

FST Antibody, HRP conjugated



PACO64240

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds directly to activin and functions as an activin antagonist. Specific inhibitor of the biosynthesis and secretion of pituitary follicle stimulating hormone (FSH).

Gene ID:

FST

Uniprot

P19883

Synonyms:

Follistatin (FS) (Activin-binding protein), FST

Immunogen:

Recombinant Human Follistatin protein (30-344AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A