

GLB1L3 Antibody



PACO64255

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:1000-1:5000

Protein Background:

vacuole, beta-galactosidase activity.

Gene ID:

GLB1L3

Uniprot

Q8NCI6

Synonyms:

Beta-galactosidase-1-like protein 3 (EC 3.2.1), GLB1L3


Immunogen:

Recombinant Human Beta-galactosidase-1-like protein 3 protein (1-236AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot. Positive WB detected in: PC-3 whole cell lysate, HeLa whole cell lysate. All lanes: GLB1L3 antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 75, 36 kDa. Observed band size: 75 kDa.