

HLA-G Antibody, FITC conjugated



PACO64269

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in the presentation of foreign antigens to the immune system. Plays a role in maternal tolerance of the fetus by mediating protection from the deleterious effects of natural killer cells, cytotoxic T-lymphocytes, macrophages and mononuclear cells.

Gene ID:

HLA-G

Uniprot

P17693

Synonyms:

HLA class I histocompatibility antigen, alpha chain G (HLA G antigen) (MHC class I antigen G), HLA-G, HLA-6.0 HLAG

Immunogen:

Recombinant Human HLA class I histocompatibility antigen, α chain G protein (25-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A