

# Hmgb1 Antibody, HRP conjugated



PACO64272

---

## Product Information

**Size:**

50ul

**Reactivity:**

Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

Hmgb1

**Uniprot**

P63159

**Synonyms:**

High mobility group protein B1 (Amphoterin) (Heparin-binding protein p30) (High mobility group protein 1) (HMG-1), Hmgb1, Hmg-1 Hmg1

**Immunogen:**

Peptide sequence from Rat Hmgb1 protein (2-18AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A